




LABOMED, INC.

www.labomed.com  
spectro@labomed.com



## LB-101 Monocular Zoom Microscope

Eyepiece:	Eyepiece WF10x
Zoom Ratio:	Zoom Ratio: 1:6.5 (0.7×-4.5×)
Total Magnification:	7×-45×
Working Distance:	90mm
Stand Focusing Range:	55mm
Stand Clips:	90mm Frosted Glass or Black/ White Plate & Clips
Illumination:	Reflect Illumination, 6V/15W Halogen Lamp with Reflector
Bottom Illumination:	6V/20W Halogen Lamp
Video System:	C Mount