




LABOMED, INC.

www.labomed.com
spectro@labomed.com



LB-112 Monocular Biological Microscope

Eyeiece:	Huygens Eyeiece H5x, H10x, H16x
Objective:	Achromatic Objective 10x, 40x, 100x
Noseiece:	Triple Noseiece
Stage:	Attachable Mechanical Stage 110-120mm / 60x30mm
Condenser:	Abbe with Iris Diaphragm and Filter
Illumination System:	Plan-concave Mirror
Filter:	Blue Filter, Green Filter