




LABOMED, INC.

www.labomed.com
spectro@labomed.com



LB-123 USB Digital Microscope

Easy to use portable digital microscope with USB output interface.

20× and 300× Magnification.

5.0 Mega Pixels CMOS sensor.

Manual Focus from 10mm to 500mm.

LED illumination with 8pcs LED lamps, brightness adjustable.

Photo management software, Auto White Balance and Exposure.

Powered by USB 2.0 cable, 5V DC.

Metal stand to hold the microscope steady.

Multi-language support: English, Germany, Spanish, French, Russian.

Operating System: Windows Win7, Win 8, Win 10 32 bit & 64 bit,
Mac OS X 10.5 or above.