




LABOMED, INC.

www.labomed.com
spectro@labomed.com



LB-282 LED Fluorescent Binocular Medical Diagnostic Microscope

Six-holes disk fluorescence unit, More fluorescence filters can be supplied

Reckless Wire Driven Mechanical Fixed Stage with focus lock is safer than traditional stage

Quintuple turret phase contrast unit with 10X/20X/40X/100X infinity plan phase contrast objective for phase contrast and bright field observation.

N.A.0.9/0.13 Swing-out Condenser

Dark field condenser (dry) available to 4X-40X Objective

Dark field condenser (wet) available to 100X Objective

Infinity Plan Objectives

This model is mainly used for sputum smear examination of Tuberculosis. The fluorescent filter for O Auramine is mainly used for this application. This microscope is mainly used in the Centers for Disease Control and Prevention (CDC).