



LABOMED, INC.

www.labomed.com
spectro@labomed.com

LB-320 Academic Zoom Stereo Microscope with LED Illumination

Features

<p>97mm</p>	<p>High resolution, Wide viewing field, Long working distance</p>
	<p>Diopter adjustable</p>
	<p>Zoom objective 0.7×-3×/ 1×-4.5×, Zoom ratio 1:4.3/ 1:4.5</p>
<p>45°</p>	<p>Binocular head inclined at 45°, Interpupillary 52-75mm</p>

Dimensions

